HOSPITAL
NACIONAL
ROSALES

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ULTRALITE

ENTERPRISES, INC.

Operator/Service Manual
For
Professional Full Body Units
With
6809 Computer Control Systems

ULTRALITE ENTERPRISES, INC.
390 FARMER COURT
LAWRENCEVILLE, GA 30046
SERVICE: 770-963-0594 Fax: 770-995-7171
ultraliteuv.com
sales@ultraliteuv.com

CONTROL LOCKED

This is the first large rectangular box in the upper left-hand corner. When the light is solid red, it informs the operator that the control system is in the locked status and cannot be operated until unlocked using the user or maintenance password. If the light is blinking, it informs the operator that the control system is in the locked status and that a key has been pressed. Before you enter your password you must press ENTER to clear the entry. The controller will display 'INVALID PASSWORD', then go back to the day and date with the red light in the solid position. At this time, you may enter the user or maintenance password in order to unlock the controller. The light will blink as the password is being entered and will turn off when ENTER/YES is pressed; assuming the correct password has been entered. The operator should always lock the controller at the end of the day or when not in use to prevent unauthorized use, or misuse.

ACTIVE CIRCUIT (Front, Rear)

These indicator lights communicate which lamp circuits, or banks, are active during a treatment. During treatment set-up, they blink until the desired circuit is selected. The active circuits will remain on; the inactive turns off. The front bank of lamps is the half of the machine that has the controller attached to it. The rear bank is the half of the machine that does not have the controller attached to it.

CALIBRATE & REFERENCE

These lights remind the operator to reference or recalibrate the unit. They are normally off indicating that the machine has been calibrated or referenced within a reasonable amount of operating hours. When these indicator lights blink, the controller is indicating that the operation hours have exceeded the recommended limit since the last calibration or reference. You should perform these operations in a timely manner or the controller will force you to do so before allowing the operator to continue. Please check pages 17, 18, 19, 20, 21, 26, 40, & 41 for more information.

DOSIMETER ACTIVE (A, B)

This light tells you that you have asked for a dose treatment expressed in Joules for either a UVA or UVB treatment when on. During setup of a treatment, the lights will blink indicating that you have the option of selecting this mode of operation. If the time treatment is selected, then the dose treatment is no longer valid and the indicator lights will turn off. If you select a dose treatment, for UVA or UVB, the timer circuit is still active, but the indicator lights will turn off. Any time after the treatment is started, the elapsed time of the treatment can be viewed using the TIMED UVA or TIMED UVB keys on the controller. To return to the dose display, press the DOSE UVA or DOSE UVB button on the controller.

TIMER ACTIVE (A, B)

This light tells you that you have asked for a TIMED treatment when illuminated. The timer functions are always active, but are secondary when dose is selected as the primary function. These lights blink during treatment set-up to indicate that you have the option of selecting the time only mode of operation. If the dose treatment is selected, then the time option is no longer valid and the light will remain off. However, any time after the treatment is started, using the dose option, the elapsed time of the treatment can be viewed using the TIMED UVA or TIMED UVB keys on the controller.

SENSOR WARNING

When this light is on, it will be preceded by a beeping sound that will indicate an error situation in the sensor circuits (during a DOSE treatment only). During treatment, there will be a sensor warning if the intensity level changes by more than 12.5% per second. This can happen if the patient gets too close to the sensor or anything blocks the view of the sensor to the lamps. If the errors persist for more than 5 seconds, the control unit will pause the treatment until the error has been corrected. When the error has been corrected, operation of the unit may be continued by pressing RESUME. The controller will continue the treatment where it stopped with no need to re-enter anything. If the beeping continues after checking everything, it is possible you may have a problem with the controller or sensor. The SENSOR WARNING light will stay on until the treatment has been completed and cleared.

HOLD (PATIENT & OPERATOR)

This light goes on when the controller or someone has paused the treatment. The patient can execute this by pressing the patient hold button located inside the unit. If the operator hold light is illuminated, then the controller has interrupted the treatment for one of three possible reasons; An error in the sensor circuit that persists for more than 5 seconds, a power interruption, or failure in some section of the controller itself. The next warning light will tell you where the hold condition originated.

SELECT CODE 63

Write non-volatile RAM. . Used only to reset entire system. **DO NOT USE** unless instructed to do so by an Ultralite service representative. This operation will clear all calibration, lamp ages, clock values, and all other memory used by the controller.

SELECT CODE 64

Displays gain values. This allows the operator to view the gain values used by the controller for dose treatments and intensity readings.

SELECT CODE 65

CHANGE PASSWORD? Enter new password. This can be used at any time to establish a new user password and to clear the previous password.

SELECT CODE 66

CONFIGURATION CODE. The configuration number indicates what operations the machine is setup to perform. Controllers are configured using the following codes:

UVA Sensor	1
UVB Sensor	2
Dosimetry	4
UVA Lamps	8
UVB Lamps	16
Printer	128

By adding the numeric values corresponding to the options installed on the machine, the controller can allow or limit the operations of different procedures. A machine having UVA Lamps, UVB Lamps, using dosimetry, and having an UVA sensor would have a configuration value of $1+2+4+8+16=29$.

SELECT CODE 67, 68, 69, 70

Allows the operator to reset the lamp ages of the four different banks of the machine. The banks are Front UVA, Rear UVA, Front UVB, and Rear UVB respectively. The lamp ages should be reset to zero each time lamps are changed in the machine. This allows the operator to keep track of the operational hours of each set of lamps.

SELECT CODE 71, 72, 73, 74, 75, 76

FACTORY USE ONLY. Used to set serial interface options and to set other modes of operation.

CALIBRATION & REFERENCING

These two words are used a lot when using Ultralite phototherapy unit and it is important that the user knows the difference between the two before we explain how to perform these functions. When performing a calibration, an independent meter will be needed to take a reading of the intensity of the lamps because selecting Code 5 will automatically erase previous stored readings. When you REFERENCE the unit you are only updating the stored intensity registers.

CALIBRATION

The Ultralite phototherapy unit is calibrated using a standard photometer, before it leaves the plant. The sensors are mounted on the interior handle and are located as such that they receive slightly more light than the patient, but they receive it in a consistent and predictable manner. The capabilities of the controller allow the system to simulate intensity levels at any reasonable location. The sensors should always be checked for any kind of obstruction of light reaching the lenses. For example, a patient leaving glasses or clothing hanging over sensor or even touching the sensor with sunscreen on their hands. It is recommended to clean the lenses with a glass cleaner regularly.

The process of calibration should be repeated occasionally to compensate for any changes or "drift" that might occur in any of the related components. This drift will be very slight and gradual. The computer has an internal section that maintains when the last calibration was performed. WHEN 500 HOURS OF USE HAS ELAPSED SINCE THE LAST CALIBRATION, THE CALIBRATE INDICATOR LIGHT WILL TURN ON AND WILL NOT TURN OFF UNTIL THE UNIT HAS BEEN CALIBRATED. This will not disable your unit. It is only a warning that it needs to be done. Ultralite keeps loaner meters to rent for this purpose and we also sell what we believe to be the finest Photometer on the market today. Please see page 18 for photometer specifications.

PH10AB PHOTOMETER SPECIFICATIONS

ULTRALITE PH10AB PHOTOMETER


SPECIFICATIONS:

IT HAS BEEN SAID, "NECESSITY IS THE MOTHER OF INVENTION". THE DEVELOPMENT OF THE PH10AB INTEGRATING PHOTOMETER IS A PERFECT ILLUSTRATION OF THIS CLICHÉ.

NON-AVAILABILITY OF A SATISFACTORY UV MEASURING DEVICE FOR PHOTOTHERAPY EQUIPMENT COMPELLED US TO DEVELOP OUR OWN. WE HAVE INTRODUCED THIS ADVANCED TECHNOLOGY TO THE PHOTOTHERAPY INDUSTRY TO ENCOURAGE STANDARDIZATION IN UV PHOTOMETRIC MEASUREMENTS. THE PH10AB CAN BE USED WITH ANY TYPE OF PHOTOTHERAPY EQUIPMENT AND INCREASES THE ABILITY FOR SAFE TREATMENTS WITH EQUIPMENT LACKING RELIABLE DOSIMETRY CONTROL.

TECHNICAL SPECIFICATIONS:

- 9 SENSORS FOR UVA AND UVB
EITHER TYPE OF UV CAN BE MEASURED
- LARGE LCD DIGITAL DISPLAY
RESULTS CAN BE EASILY READ
- PORTABLE
- COMPACT AND SELF-CONTAINED
- EVERYTHING NEEDED IS STORED IN A DURABLE CASE
- MICROPROCESSOR CONTROLS
STATE-OF-THE ART TECHNOLOGY IS USED FOR CONTROL
- MEMORY FOR STORING READINGS
ALL THE DATE FROM READINGS CAN BE STRED FOR RECALL LATER
- FLOATING DECIMAL
YOU CAN CHOOSE YOUR DECIMAL PLACEMENT
- TWO RESPONSE SPEEDS
INSTANT READ-OUTS OR AVERAGE VALUES CAN BE SELECTED
- INTRGRATING DOSIMETRY
NO CHARTS OR CALCULATIONS ARE NEEDED FOR FIGURING DOSE OR INTENSITY

WARRANTY: ONE-YEAR FULL WARRANTY WITH ONE CALIBRATION INCLUDED

DIMENSIONS:

9.5" X 11.5" X 2.75"

7 LBS WEIGHT APPROXIMATE

STEPS TO CALIBRATE:

1. Turn the meter on with the "On" key (red button on side). The display will show the software version number for one second in the Intensity display. Then the serial number will appear in the Intensity display, and the last calibration date in the Dose display. If the meter was last calibrated more than one year ago, recalibration is recommended before use.
2. With the sensors held in the photometer case (shielding them from all light input), select the mode (UVA or UVB). If the Intensity shows more than 0.05 mW/cm^2 , the meter should be recalibrated before use.
3. Clip the sensors to the provided mount and use the UVA or UVB key to select the appropriate display.
4. Be sure that the meter is in the "Slow" mode (indicator on in the "Slow" window). If the indicator is off, use the "Speed" key to activate it.
5. The intensity display shows the average photometric intensity during the last 10 seconds.
6. Unlock the 6809 controller by first pressing the red space/lock button then type Ultramain, and press ENTER.
7. A. Press SELECT, 5, ENTER. The display will read "CALIBRATION REQD" for a few seconds, and then display "SENSORS COVERED" At this time, just be sure no ultraviolet lights are on in the room and the door of the machine is closed.
B. If "FUNC INT CALC" is displayed, instead of "CALIBRATION REQD", Then press "SELECT, 5, ENTER. Configuration will then be displayed. Put in your machine configuration #, then press ENTER.
C. Re-enter SELECT, 5, ENTER, CALABRATION REQD, should now be displayed.
8. Press "YES" in answer to the query "SENSORS COVERED". The UVA lamps will immediately turn on and the display will read "A SENSOR READY?" Observe the output on the meter's scale and notice that the output is unstable as the lamps are still warming to their optimum operating temperature. It may take approximately 5 or 6 minutes.
9. When the reading has become stable and stops climbing on the scale press "ENTER". The display will now read "UVA FLUX": it is asking you for the new reading. Enter the value that your meter reads and press "ENTER". Example: To enter 29.5 mw/cm^2 , press 2, 9.5, enter.
10. The display will read "GAIN=*****29.5" for a few seconds and go directly into calibration. If unit is not equipped with a B sensor, the display will read "FUNC INT CALC" for a few seconds and turn on the UVA lamps. Press "ENTER" and the calibration is complete with the 29.5 being your new intensity level. IF YOUR UNIT IS EQUIPPED WITH UVB SENSOR, CHANGE THE METER TO ITS UVB SENSOR AND SWITCH THE DIAL TO THE UVB RANGE. MOUNT THE SENSOR IN UNIT. NOW YOU ARE READY TO CALIBRATE THE UVB.
11. If your unit is equipped with a UVB sensor it will go directly into "B SENSOR READY" as soon as you enter the intensity for UVA. Press the UVB button on the photometer. Observe the output on the meter's scale and notice that the output is unstable as the lamps are still warming to their optimum operating temperature. It may take approximately 3 minutes.
12. When the reading has become stable and stops climbing on the scale press "Enter". The display will now read "UVB Flux". Enter the value that your meter reads and press "Enter". A gain value will display momentarily and then go back into showing the intensity for UVA. Press "Enter" for both UVA and UVB.
13. **BE SURE AND RECORD YOUR NEW INTENSITY LEVELS!** All intensity values used with this computer are expressed in milli-watts per square centimeter for both UVA and UVB. Press select, 61, to read the stored intensity values. Then press clear, followed by select 64 to read and record the corrected gain values.

REFERENCING

It is important that you reference or update the data stored in the two functional intensity registers often and records those outputs. This data is related to the lamps and their on-time and not any other components of the machine. Since lamps have only a few thousand hours of useful life and degrade with use, data relating to these lamps should be updated frequently. It is done without a photometer and is a nearly automatic process of the machine.

STEPS to REFERENCE (select code 4)

Is important that you reference the data stored in the memory of the two functional intensity registers often and keep records of the outputs. Once a month is good. Once a week is better. The computer monitors the elapsed on-time and activates the "reference" indicator light after 10 hours of use. This is not a good thing. Do not let the computer go that long without referencing because it will cause the computer to automatically turn on and start referencing on it's own one day when you unlock the controller. It is much better to reference it when you have the time to let it run and warm up properly and record all the numbers. The process takes about 5-10 minutes.

1. Unlock controller with ULTRAMAIN code.
2. Press, "SELECT, 4, ENTER". The control will then turn on the UVA lamps and display will read "UVA FLUX=xx.xx". Allow lamps to reach a peak, and when numbers become stable. 5-6 minutes usually. The last three numbers will start to go up and down. Do not let it run too long or numbers will continue to fall.
3. Press "ENTER". Record the number that was on the display when you pressed ENTER. This is your new flux intensity value. Note, if you only have UVA lamps in your unit, then you are done.
4. However, if you have a combination unit with UVB lamps and no B sensor the UVB lamps will turn on but there will be no numbers. You must enter a value for the UVB intensity. Use the number 1. Press 1 and ENTER. All lamps will then turn off and the controller will go back to time and date.
5. Should you have a combination unit with a UVB Sensor attached, the UVB lamps will turn on and the display will show the UVB intensity value. Let these lamps warm up until the value becomes stable, then press ENTER.

We recommend that this is done on a regular basis, for the good of your patients and your machine. Be sure to keep a log of the reference date, the new flux value (Sheet 61), and the present corrected gain values (Select 64). Note, this will be valuable information should your machine ever lose its memory or calibration.

COMPUTATION OF CALIBRATION VALUES


Precise measurement of light intensities in phototherapy devices is extremely difficult to perform, and in fact is impossible to achieve with equipment available to most equipment users. However, if measurement and calibration are performed in a consistent and logical manner meaningful measurements can be derived and safe and effective treatment of patients can be maintained.

UVA Calibration is pretty straightforward because it is a simple matter to place a UVA detector in a position in the chamber that replicates light intensity at the skin surface. For the most precise calibration the sensor should be placed on the skin surface of a person in the machine at about chest height. This method is not always so practical either, because the clothes worn by the person holding the sensor and even the size of the person all affect the reading substantially. We have found that in Ultralite machines mounting the sensor on the center grab handle facing the rear of the machine near the hinge provides a consistent measurement that is approximately equal to intensity levels measured at the chest skin surface of most patients.

UVB Calibration is more difficult. The first consideration is that the true light intensities of UVB in a typical chamber are 1/20 to 1/40 of the UVA intensities. With such low light levels it is not possible to place a sensor in a location to detect light from a wide area as is done with UVA. Therefore it is necessary to select one lamp as a representative of all the UVB lamps in machine. For this reason it is important that all UVB lamps in a machine be of the same age and from the same manufactured batch. It is also important that all measurements be made at the same distance from the lamp. Usually about 6 inches (+/-150mm) from sensor to lamp surface is used.

The second consideration is that light sensors are very directional and detect light from a relatively small conical zone. Therefore, because UVB lamps are spread so far apart the sensor attached to independent photometers can detect energy from only some lamps at a time. When calibrating a light machine each UVB lamp should first be checked with the independent meter to insure that all lamps are approximately equal in performance. The strongest lamp should be located in the sensor's field of vision.

The third consideration is the actual number of UVB lamps in the machine. Experimentation has shown that the method of calibration currently used for photometers yields a reasonably accurate measurement of UVB intensity in Ultralite machines with 8 UVB lamps. In machines with more than eight lamps the measurement must be factored by the increased number of lamps. This is because the sensor can "see" only some of the lamps. A machine with eight lamps will appear to the photometer to be equal to one with twelve or sixteen lamps because of this limited viewing angle. It is important that the reading from the photometer is multiplied by the number of UVB lamps in the machine divided by eight. Thus, if a machine has 12 lamps, the true reading will be 1.5 times the meter reading. And with 16 lamps the true reading will be double the meter reading.

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